TA-3Dprint model

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

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Produkt/ Produkttyp

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Instructions for use

Product/product type

TPO-free, light-curing acrylic-based resin for the production of customized, 3D-printed dental models. Suitable for the production of stump models, models for crown and bridge technology, veneers, inlay and onlay models, model cast splints and aligner models, as well as for the production of implant models, mock-up and wax-up models, and situation models. TA-3Dprint model can be used in all 3D printers with DLP technology with LED 405 nm or UV-LED 385 nm.

Material type

Light-curing acrylic-based resin

Users

Professional users

Target group

Dental technicians, dentists, dental professionals

Indication/intended use

TPO-free, light-curing acrylic-based resin for producing dental models of all geometries. TA-3Dprint model can be used in all 3D printers with DLP technology with LED 405 nm or UV LED 385 nm.

Material properties/technical data

TPO-free, light-curing acrylic-based plastic for 3D printing using the DLP process with LED 405 nm or with UV LED 385 nm

Processing

Mixing


Shake the product in its original packaging for at least five minutes. If the material is already in the material tray, stir thoroughly with a spatula.

Built up process


Create the design of the construction and the support structure using CAD software, the print jobs in compliance with the machine and material parameters. Then carry out the printing process in accordance with the printer's operating instructions.

Removing the object from the build platform

Use a spatula to remove the finished printed object and any polymerized residue from the build platform.

Cleaning the printed object

The printed objects can be cleaned in an automatic cleaning unit by selecting the cleaning parameters appropriate for the material. The manufacturer's instructions for use must be observed. Carefully blow off any remaining cleaning fluid with compressed air.


TA-3Dprint model

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

Cleaning alternative with ultrasonic device

step 1:

Clean the printed objects in a separate container with isopropanol for a maximum of three minutes in an ultrasonic bath.
Carefully brush openings, holes, and crevices with a soft brush.

step 2:

Clean the printed objects in a separate container with fresh isopropanol for a maximum of three minutes in an ultrasonic bath (caution: the total cleaning time must not exceed 10 minutes!). To dry, carefully blow off the printed objects with compressed air.

Alternative:

Butyldiglycol can be used as an isopropanol-free alternative for cleaning in an ultrasonic bath. After use, the printed object must be rinsed with clean water. Please follow the manufacturer's instructions for use.

Aushärtung

Cure the cleaned print objects in a light-curing device using the post-exposure parameters appropriate for the material. (Follow the manufacturer's instructions for use.)

post-light curing alternative

Post-light curing of the printed objects in a xenon light curing unit with suitable post-light curing parameters appropriate for the material (follow the manufacturer's instructions for use).

completion

Post-process the finished printed objects using conventional methods with dental rotating tools.

Note

The product is designed for use in the dental field and must be used in accordance with the instructions for use. The manufacturer is not liable for damage resulting from any other use. In addition, it is the responsibility of the user to check the material for suitability and usability for the intended purpose, especially if this is not mentioned in the instructions for use.

Safety instructions

Liquid and vapor are highly flammable.

Inhalation:

Avoid inhaling dust, smoke, gas, mist, vapor, and aerosol, and ensure good ventilation.

Protective clothing:

Wear protective gloves, protective clothing, eye protection, and face protection, and avoid contact with uncured material. Remove contaminated clothing immediately and wash before reuse.

Skin contact:

If the material comes into contact with the skin, wash the affected area with plenty of water. If a rash or irritation occurs, seek medical advice immediately.

TA-3Dprint model

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

Eye contact:

After eye contact, rinse the affected eye for several minutes under running water with the eyelid open. If eye irritation persists, seek medical advice immediately.

Storage

The product must be stored in a dry place, protected from light and in its original packaging at room temperature (between 4°C and 28°C). Even minimal exposure to light can trigger polymerization. The packaging must be carefully resealed after use. Once the expiration date has passed, the material properties and handling are no longer guaranteed. See the label for the shelf life.

Disposal

The unpolymerized material is harmful to the environment and must therefore be prevented from entering the environment. Polymerized material is not harmful to the environment.

Dispose of residual waste in liquid form and packaging in accordance with local/regional/national/international regulations.

Explanation of the labels on the packaging

catalog number/item number

batch / lot number

observe the safety instructions in the operating instructions and user manual.

www.wp-ifu.de/3dliquid

path to the electronic user manual

paper user manual available within 7 days

manufacturer

best before

temperature range

wave length

store away from sunlight